



## FIRST PROTOTYPE INSTALLED



We are proud to announce that we with support from Kapacitet A/S this weekend successfully installed our first fully functional prototype - the P3 system – at Aarhus University Hospital in Denmark.

The P3 system is capable to convert  $^{15}\text{O}$ -water gas from the cyclotron into  $^{15}\text{O}$ -water and inject it into the patient undergoing a PET exam.

The P3 system will now undergo clinical validation and animal testing before being made available for research.



### Further Information:

Martin Stenfeldt, CEO MedTrace Pharma A/S  
+45 28 10 41 49 (GMT+1)  
+1 (833) 615-4261 (USA)  
martin@medtrace.dk